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Sample	Jun(MeanCt)	GAPDH(MeanCt)	E
处理前	18	17	0.95
处理后	16	17.4	0.95

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TIANGEN BIOTECH (BEIJING) CO.,LTD

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The Plant Cell, Vol. 21: 1495–1511, May 2009, www.plantcell.org © 2009 American Society of Plant Biologists

***Arabidopsis* ASA1 Is Important for Jasmonate-Mediated Regulation of Auxin Biosynthesis and Transport during Lateral Root Formation**

Jiaqiang Sun,^{a,b} Yingxiu Xu,^{a,c} Songqing Ye,^d Hongling Jiang,^a Qian Chen,^{a,b,c} Fang Liu,^{a,b,c} Wenkun Zhou,^{a,c} Rong Chen,^{a,b,c} Xugang Li,^{b,e} Olaf Tietz,^{b,e} Xiaoyan Wu,^a Jerry D. Cohen,^d Klaus Palme,^{b,e} and Chuanyou Li^{a,b,1}

^a State Key Laboratory of Plant Genomics, National Centre for Plant Gene Research, Institute of Genetics and Developmental Biology, Chinese Academy of Sciences, Beijing 100101, China

^b Chinese-German Joint Group for Plant Hormone Research, Institute of Genetics and Developmental Biology, Chinese Academy of Sciences, Beijing 100101, China

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Journal of Experimental Botany, Page 1 of 12
doi:10.1093/jxb/erp223

RESEARCH PAPER

Molecular mechanism for jasmonate-induction of anthocyanin accumulation in *Arabidopsis*

Xiaoyi Shan^{1,*}, Yongsheng Zhang^{1,*}, Wen Peng¹, Zhilong Wang² and Daoxin Xie^{1,†}

¹ MOE Key Centre of Bioinformatics, Department of Biological Sciences and Biotechnology, Tsinghua University, Beijing 100084, China

² Temasek Life Science Laboratory, National University of Singapore, 117604, Republic of Singapore

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




Minimum Information for Publication of Quantitative Real-Time

PCR Experiments <http://miqe.info/>

Stephen Bustin BA(Mod) PhD (TCD) FSB
Professor of Molecular Science
Queen Mary University of London

MIQE

[Home](#)[News](#)[Data Standard](#)[Development](#)[MIQE](#)[Tools](#)[Consortium](#)


Real-time PCR Data Markup Language

Tools

Create RDML file
Validate RDML file
Browse libraries

Join or support

The Real-time PCR Data Markup Language (RDML) is a structured and universal data standard for exchanging quantitative PCR (qPCR) data. The data standard should contain sufficient information to understand the experimental setup, re-analyse the data and interpret the results. The data standard is a compressed text file in Extensible Markup Language (XML) and enables transparent exchange of annotated qPCR data between instrument software and third-party data analysis packages, between colleagues and collaborators, and between authors, peer reviewers, journals and readers. To support the public acceptance of this standard, both an on-line RDML file generator is available for end users, as well as RDML software libraries to be used by software developers, enabling import and export of RDML data files.

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The Minimum Information for Biological and Biomedical Investigations

-  Portal
-  Foundry
-  MICheckout
-  MIBBI Search
-  MIBBI News
-  Discussion
-  Related
-  About

Project News Highlights

- The report on the last [MIBBI workshop](#) is [available free](#) from the journal [Standards in Genomic Sciences](#).
- MIBBI has a new **funding stream** as part of the [BioSharing](#) network, supported by BBSRC and NERC funds to the [ISA project](#).
- MIBBI Portal use recommended by **BMC journals** ([example](#)) and **OMICS** ([editorial](#))
- **Free download:** [The MIBBI paper \(Nature Biotechnology\)](#) & [supplementary information \(additional figures\)](#)

Site Navigation

[The MIBBI Portal](#)

Access to Minimum Information guidelines from diverse bioscience communities

[MICheckout](#)

Compile made-to-measure checklists in a variety of formats, or just for display.

[The MIBBI Foundry](#)

Towards the next generation of MI guidelines for the biosciences

[MIBBI search](#)

A Google™ Custom Search Engine covering a range of relevant web sites.

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Clinical Chemistry 55:4
000–000 (2009)

Reviews

The MIQE Guidelines: *Minimum Information for Publication of Quantitative Real-Time PCR Experiments*

Stephen A. Bustin,^{1*} Vladimir Benes,² Jeremy A. Garson,^{3,4} Jan Helleman,⁵ Jim Huggett,⁶
Mikael Kubista,^{7,8} Reinhold Mueller,⁹ Tania Nolan,¹⁰ Michael W. Pfaffl,¹¹ Gregory L. Shipley,¹²
Jo Vandesompele,⁵ and Carl T. Wittwer^{13,14}

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ITEM TO CHECK	IMPORTANCE	CHECKLIST
EXPERIMENTAL DESIGN		
Definition of experimental and control groups	E	
Number within each group	E	
Assay carried out by core lab or investigator's lab?	D	
Acknowledgement of authors' contributions	D	
SAMPLE		
Description	E	
Volume/mass of sample processed	D	
Microdissection or macrodissection	E	
Processing procedure	E	
If frozen - how and how quickly?	E	
If fixed - with what, how quickly?	E	
Sample storage conditions and duration (especially for FFPE samples)	E	
NUCLEIC ACID EXTRACTION		
Procedure and/or instrumentation	E	
Name of kit and details of any modifications	E	
Source of additional reagents used	D	
Details of DNase or RNase treatment	E	
Contamination assessment (DNA or RNA)	E	
Nucleic acid quantification	E	
Instrument and method	E	
Purity (A260/A280)	D	
Yield	D	
RNA integrity method/instrument	E	

Table 1. MIQE checklist for authors, reviewers, and editors.^a

Item to check	Importance	Item to check	Importance
Experimental design		qPCR oligonucleotides	
Definition of experimental and control groups	E	Primer sequences	E
Number within each group	E	RTPrimerDB identification number	D
Assay carried out by the core or investigator's laboratory?	D	Probe sequences	D ^d
Acknowledgment of authors' contributions	D	Location and identity of any modifications	E
Sample		Manufacturer of oligonucleotides	D
Description	E	Purification method	D
Volume/mass of sample processed	D	qPCR protocol	
Microdissection or macrodissection	E	Complete reaction conditions	E
Processing procedure	E	Reaction volume and amount of cDNA/DNA	E
If frozen, how and how quickly?	E	Primer, (probe), Mg ²⁺ , and dNTP concentrations	E
If fixed, with what and how quickly?	E	Polymerase identity and concentration	E
Sample storage conditions and duration (especially for FFPE ^b samples)	E	Buffer/kit identity and manufacturer	E
Nucleic acid extraction		Exact chemical composition of the buffer	D
Procedure and/or instrumentation	E	Additives (SYBR Green I, DMSO, and so forth)	E
Name of kit and details of any modifications	E	Manufacturer of plates/tubes and catalog number	D
Source of additional reagents used	D	Complete thermocycling parameters	E
Details of DNase or RNase treatment	E	Reaction setup (manual/robotic)	D
Contamination assessment (DNA or RNA)	E	Manufacturer of qPCR instrument	E
Nucleic acid quantification	E	qPCR validation	
Instrument and method	E	Evidence of optimization (from gradients)	D
Purity (A_{260}/A_{280})	D	Specificity (gel, sequence, melt, or digest)	E
Yield	D	For SYBR Green I, C_q of the NTC	E
RNA integrity: method/instrument	E	Calibration curves with slope and y intercept	E
RIN/RQI or C_q of 3' and 5' transcripts	E	PCR efficiency calculated from slope	E
Electrophoresis traces	D	CI _s for PCR efficiency or SE	D
Inhibition testing (C_q dilutions, spike, or other)	E	r^2 of calibration curve	E
Reverse transcription		Linear dynamic range	E
Complete reaction conditions	E	C_q variation at LOD	E
Amount of RNA and reaction volume	E	CI _s throughout range	D
Priming oligonucleotide (if using GSP) and concentration	E	Evidence for LOD	E
Reverse transcriptase and concentration	E	If multiplex, efficiency and LOD of each assay	E
Temperature and time	E	Data analysis	
Manufacturer of reagents and catalogue numbers	D	qPCR analysis program (source, version)	E

C _q s with and without reverse transcription	D ^c	Method of C _q determination	E
Storage conditions of cDNA	D	Outlier identification and disposition	E
qPCR target information		Results for NTCs	E
Gene symbol	E	Justification of number and choice of reference genes	E
Sequence accession number	E	Description of normalization method	E
Location of amplicon	D	Number and concordance of biological replicates	D
Amplicon length	E	Number and stage (reverse transcription or qPCR) of technical replicates	E
In silico specificity screen (BLAST, and so on)	E	Repeatability (intraassay variation)	E
Pseudogenes, retropseudogenes, or other homologs?	D	Reproducibility (interassay variation, CV)	D
Sequence alignment	D	Power analysis	D
Secondary structure analysis of amplicon	D	Statistical methods for results significance	E
Location of each primer by exon or intron (if applicable)	E	Software (source, version)	E
What splice variants are targeted?	E	C _q or raw data submission with RDML	D

^a All essential information (E) must be submitted with the manuscript. Desirable information (D) should be submitted if available. If primers are from RTPrimerDB, information on qPCR target, oligonucleotides, protocols, and validation is available from that source.

^b FFPE, formalin-fixed, paraffin-embedded; RIN, RNA integrity number; RQI, RNA quality indicator; GSP, gene-specific priming; dNTP, deoxynucleoside triphosphate.

^c Assessing the absence of DNA with a no-reverse transcription assay is essential when first extracting RNA. Once the sample has been validated as rDNA free, inclusion of a no-reverse transcription control is desirable but no longer essential.

^d Disclosure of the probe sequence is highly desirable and strongly encouraged; however, because not all vendors of commercial predesigned assays provide this information, it cannot be an essential requirement. Use of such assays is discouraged.

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Absolute RNA® Purification Kits
Absolute RNA® FFPE Kit (石蜡包埋样品)
miRACLE miRNA Isolation Kit (miRNA提取)
SideStep® QPCR MasterMix Kit (无需RNA提取)

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QPCR 扩增

Brilliant® II QPCR and QRT-PCR Reagent
High-specificity miRNA QRT-PCR Detection Kit

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

APP for cell-phone

Real-time PCR and Next Generation Sequencing
Changing the world of molecular diagnostics
and molecular research

November 13th-16th, 2012 - Daejeon, South Korea
SYMPOSIUM - EXHIBITION - WORKSHOPS

[Home](#)

NEWS

- Changed dates!
The qPCR symposium 2012 will be November 13th-16th
- The qPCR symposium 2012 will be November 6th-9th in Daejeon, South Korea

TATAA Biocenter, PhileKorea Technology and Select Biosciences proudly presents and welcomes you to the 1st International Genomics Research Asia

The conference will consist of two days with scientific sessions and industrial exhibition and also conference workshops. Besides interesting talks about the new developments in real-time qPCR and Next Generation Sequencing (NGS) research, special focus in this conference is:

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


The interface includes several panels for parameter adjustment and analysis:

- Sequence Parameters:** Includes fields for File Name, Length (1310 bp), and a sequence viewer.
- Primer Parameters:** Includes checkboxes for Primer Tm, Primer GC Content, Primer Length, Primer Composition, and Primer Secondary Structure, each with associated sub-parameters.
- Primer Selection:** A list of candidate primers with checkboxes for selection.
- Primer Design Options:** Radio buttons for Forward Primer, Reverse Primer, and Probe.
- Warning and Comments:** A large text area for displaying warnings and user comments.
- Export Options:** Buttons for Export, All Store, and Export Order Information.

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

OLIGO 6

File Edit Analyze Search Select Change View Window Help

☐ Display Hairpin with 2 or more bp Stems

Upper Primer Mus beta actin:954U20

The most stable 3'-dimer: 3 bp, -3.2 kcal/mol

5' ACCTCTATGCCAACACAGTG 3'

3' GTGACACAACCGTA

The most stable dimer overall: 3 bp, -3.2 kcal/mol

5' ACCTCTATGCCAACACAGTG 3'

3' GTGACACAACCGTA

Optimal Annealing Temperature: 52.9° (Max: 67.2°)

	Position and Length	T _m [°C]	GC [%]	P.E.#
Product	92	83.8	53.3	-----
Upper Primer	954 20	64.7	50.0	375/375
Lower Primer	1026 20	64.2	50.0	377/377

Product T_m - Lower Primer T_m: 19.5

Primers T_m difference: 0.4

Concentration

Upper Primer: 200.0 nM

Lower Primer: 200.0 nM

Monovalent Cation: 50.0

Free Mg [2+]: 0.7

Total Na[+] Equivalent: 155.8

☐ Display Hairpin with 2 or more bp Stems

Lower Primer Mus beta actin:1026L20

The most stable 3'-dimer: 3 bp, -3.2 kcal/mol

5' AGGAGCCAGAGCAGTAATCT 3'

3' TCTAATGACGAGACCGAGGA

The most stable dimer overall: 3 bp, -3.2 kcal/mol

5' AGGAGCCAGAGCAGTAATCT 3'

3' TCTAATGACGAGACCGAGGA

Mus beta actin - U...

Upper Primer Mus beta actin:954U20

T_d = 59.9° [nearest neighbor method]

T_m = 73.2° [%GC method]

T_m = 60° [2(A+T)° + 4(G+C)° method]

T_m(RNA)[1M Na] = 84.2° [%GC method]

T_m(DNA:RNA)[1M Na] = 78.2° [%GC method]

A₂₆₀/A₂₈₀ = 1.88 [single strand]

M_r = 6.1K [one strand]

M_r = 12.4K [two strands]

μ g/OD = 47.4 [dsDNA]

Base	Number & %
A	6 [30.0%]
C	7 [35.0%]
G	3 [15.0%]
T	4 [20.0%]
A + T	10 [50.0%]
G + C	10 [50.0%]

Mus beta actin - Upper...

☒ Display Hairpin with 2 or more bp Stems

Primers: Mus beta actin:954U20 Mus beta actin:1026L20

The most stable 3'-dimer: 2 bp, -1.9 kcal/mol

5' ACCTCTATGCCAACACAGTG 3'

3' TCTAATGACGAGACCGAGGA

The most stable 3'-dimer: 2 bp, -1.6 kcal/mol

5' ACCTCTATGCCAACACAGTG 3'

3' TCTAATGACGAGACCGAGGA

The most stable dimer overall: 3 bp, -5.0 kcal/mol

5' ACCTCTATGCCAACACAGTG 3'

3' TCTAATGACGAGACCGAGGA 5'

Mus beta actin - L...

Lower Primer Mus beta actin:1026L20

T_d = 60.4° [nearest neighbor method]

T_m = 73.2° [%GC method]

T_m = 60° [2(A+T)° + 4(G+C)° method]

T_m(RNA)[1M Na] = 84.2° [%GC method]

T_m(DNA:RNA)[1M Na] = 78.2° [%GC method]

A₂₆₀/A₂₈₀ = 2.04 [single strand]

M_r = 6.3K [one strand]

M_r = 12.4K [two strands]

μ g/OD = 47.4 [dsDNA]

Base	Number & %
A	7 [35.0%]
C	4 [20.0%]
G	6 [30.0%]
T	3 [15.0%]
A + T	10 [50.0%]
G + C	10 [50.0%]

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OLIGO - Mus beta actin

File Edit Analyze Search Select Change View Window Help

for Primers and Probes F3
for Hairpin Loop Stems
for a Sequence String
for Palindromes
for Restriction Sites...
for Restriction Sites in Protein...

Memory Tables
Copy Memory Table...

Temperature: 52.9° (Max: 67.2°)

Position and Length Tm [°C] GC

92 83.8 53

50

50

Mus beta actin - Primer P...

[1] [2] Sort Sort Sort

#	Positions of Primers	Prod. Len.	Opt. Ta	%GC
1	43 244	221	60.4	63.8
2	43 607	584	60.3	60.1
3	44 243	219	60.0	64.4
4	44 244	220	61.0	64.1
5	44 245	221	59.9	64.3
6	44 245	221	61.0	64.3
7	44 412	388	59.6	60.8
8	44 607	583	60.9	60.2
9	44 832	808	60.7	59.0
10	44 833	809	60.7	59.1
11	44 1404	1380	58.0	52.9
12	45 243	218	59.9	64.2
13	45 409	384	58.5	61.2
14	45 410	385	58.4	61.0
15	45 411	386	58.4	60.9
16	45 412	387	59.6	60.7
17	45 1202	1177	57.4	56.9
18	45 1203	1178	57.5	56.9
19	45 1402	1377	56.7	52.8
20	45 1403	1378	58.0	52.8
21	45 1403	1378	56.7	52.8
22	45 1404	1379	58.0	52.9
23	46 409	383	59.4	61.1
24	46 410	384	59.4	60.9
25	46 411	385	59.3	60.8
26	46 826	800	59.3	58.9
27	46 1202	1176	57.4	56.9

Search for Primers and Probes

☒ + Strand Search ☒ - Strand Search

PCR Primers: ☐ Consensus Primers ☒ Complex Substrate

☒ Compatible Pairs
☐ Compatible with the Upper Primer
☐ Compatible with the Lower Primer

☐ Sequencing Primers
☐ Hybridization Probes

☒ Duplex-free Oligonucleotides
☒ Highly Specific Oligos (3'-end Stability)
☒ Oligonucleotide with GC Clamp
☒ Eliminate False Priming Oligonucleotides
☐ and Continue Above Search in Other File(s)
☒ Oligonucleotides within Selected Stability Limits
☒ Hairpin-free Oligonucleotides
☒ Eliminate Homooligomers / Sequence Repeats
☒ Eliminate Frequent Oligos

Search Mode:
☒ Mark
☐ Unmark

OK Cancel Search Ranges Parameters Defaults

Base Number & %

Base	Number	%
A	6	[30.0%]
C	7	[35.0%]
G	3	[15.0%]
T	4	[20.0%]
A + T	10	[50.0%]
G + C	10	[50.0%]

Search for sequencing and PCR primers or hybridization probes

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发表文章	杂志
Functional screen reveals SARS coronavirus nonstructural protein nsp14 as a novel cap N7 methyltransferase	PNAS
Arabidopsis ASA1 Is Important for Jasmonate-Mediated Regulation of Auxin Biosynthesis and Transport during Lateral Root Formation	PLANT CELL
The Arabidopsis RING Finger E3 Ligase RHA2a Is a Novel Positive Regulator of Absciscic Acid Signaling during Seed Germination and Early Seedling Development	Plant Physiology
Molecular mechanism for jasmonate-induction of anthocyanin accumulation in Arabidopsis	J. Exp. Bot
Glyceraldehyde-3-phosphate dehydrogenase of Xanthomonas campestris pv. campestris is required for extracellular polysaccharide production and full virulence	Microbiology
Toll-like receptor 2 mediates the induction of IL-10 in corneal fibroblasts in response to Fusarium solu	Immunology and Cell Biology
Inhibition of telomerase RNA (hTR) in cervical cancer by adenovirus-delivered siRNA	Cancer Gene Therapy
Local hyperthermia induces apoptosis of keratinocytes in both normal skin and condyloma acuminata via different pathways	Apoptosis
Spinal D-Amino Acid Oxidase Contributes to Neuropathic Pain in Rats	J. Pharmacol. Exp. Ther
Smad2/3 activities are required for induction and patterning of the neuroectoderm in zebrafish	Developmental Biology
Cloning the PvP5CS gene from common bean (Phaseolus vulgaris) and its expression patterns under abiotic stresses	Journal of Plant physiology

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